




sartorius stedim
biotech

Discover our convenience packages

no operational costs for 2 years with the UniVessel® Consumable Kit*



* depending on frequency of use

Fermentor|Bioreactor + Consumable Kit

Get a consumable kit for free



BIOSTAT® Aplus

The BIOSTAT® Aplus is a compact, autoclavable fermentor|bioreactor system specially designed for educational use and preliminary or investigational R&D applications. The application-driven, configured packages for microbial and cell culture include everything needed for to get started immediately.



BIOSTAT® Bplus

The BIOSTAT® Bplus is a compact autoclavable laboratory fermentor | bioreactor with graphical user interface and touch screen. The integrated system concept enclose in a stainless steel housing the control-, thermostat-aeration system, a stirrer control unit as well as up to four peristaltic pumps. The system satisfy the demands of both microbial & cell culture applications.



Free Consumable Kit

The Consumable Kit helps you to ensure the operability of your fermentor|bioreactor for around 2 years (depending on frequency of use)

Content of the Consumable Kit

- Spare glass vessel
- Spare pH probe
- DO membrane kit
- Consumable Kit UniVessel®
- Mediafilter Sartolab P20plus
- Midisart filter

Fax Reply Form
+49.5661.929949

Name	
Company	
Department	Position
Address	
Post or Zip Code	City Country
Phone	Fax
E-Mail	

Sartorius Stedim Systems GmbH
Phone +49.5661.71.3400
info.fermentation@sartorius-stedim.com
www.sartorius-stedim.com

Promotion valid until 31th March, 2010

Promotion Code: Consumable Kit